



Key Features

- Audible and Visual Alarms
- Three Alarm Functions
- Air Cooling,Auto defrost

.....

More Detail:

Specifications

Technical Data

Type	Under-counter Pharmacy Refrigerator
Cabinet Type	Chest
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise Lever (dB)	30

Performance

Temperature Range(°C)	2~8
------------------------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
Power(W)	160
Electrical Current(A)	1.05

Dimensions

Capacity(L/Cu.Ft)	118/4.2
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	512*415*638(mm) 20.2*16.3*25.1(in)

Exterior Dimensions(W*D*H)	597*645*835(mm) 23.5*25.4*32.9(in)
Packing Dimensions(W*D*H)	675*680*930(mm) 26.6*26.8*36.6(in)
Container load (20'/40'/40'H)	48/102/102

Functions

High/Low Temperature	Y
Remote Alarm	N/A
Power Failure	N/A
Sensor Error	Y
Low Battery	N/A
Door Ajar	Y

Accessories

Caster	N/A
Foot	Y
Test Hole	N/A
Shelves/Inner Doors	3/-

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate	CE
--------------------	----